GIMA

TUBI SILICONE
SILICONE TUBE
TUYAU EN SILICONE
TUBO DE SILICONA
TUBO DE SILICONE
ΣΩΛΗΝΑΣ ΣΙΛΙΚΟΝΗΣ

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ

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SILICONE TUBE

The silicone tube has a wide range of use in the medical field, its great stability at high and low temperature (-50/+250°C) makes it suitable for several uses. Its typical employment is the liquids or gas conduction, but only and exclusively for EXTERNAL use.

Technical specifications

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Colour: transparent

Strength: 60° ±5% Shore - DIN 53505

Density: 1.17 g/ccm - DIN 53479A

Tension: 10.8 N/mm² - DIN 53504 S1

Extension: 490/530% - DIN 53504 S1

Traction resistance: 34-35 N/mm - ASTM D 624 B

Resilience: 51% - DIN 53512

Compression set: (22h/175°C): 30-15% - DIN 53517

Infrared rays vulcanized

Treated at a temperature of 200°C for 10 hours

Non-sticking, it avoids encrustments and coagulations.

- In addition to the notes in the technical specifications, please check the compatibility of the used fluids and any possible deterioration in the course of time.
- Check that no sharp edge or corner that can cause injuries have been produced, while cutting the tube.
- Check the tube tightness, particularly after several reuses.
- The dimensions of the tube (inner diameter, external diameter and length) must fit the connections and/or the opening where the tube has to be inserted.

- *Do not use in a not proper way*
- *Do not arrange it in an unsuitable way, so to cause danger or impediment.*
- *Do not knot to close it at one end*
- *The item is sold NON sterile. Sterilize before using!*

The silicone tube can be reused and sterilized. The sterilization of the item can be carried out: in autoclave (135° C – 2.2-2.5 bar – please see details about autoclave), through gamma rays or hot steam at a temperature of 200°C.

The disposal of the item must meet the local standards concerning dangerous medicines waste.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product.

This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA.