

FILTRI BATTERIOLOGICI SPIROGUARD conf. 50 pz.



FILTRO ANTIBATTERICO SPIROGUARD - conf. da 50 pz.

Compatibilità di Spiroguard:

Erich Jaeger Gmbh Masterlab e Masterscreen

Ferraris Medical - P.K. Morgan Autolink-Benchmark Transflo e USA Transfer System

Sensormedics 2200 - 2400 - 2450

Vitalograph Wedge-Bellow-Alpha e Escort Spirometers