

**CEROTTI SENSITIVE 5 misure - 80 conf. da 40 pz.
conf. 3200 pz.**



Cerotti ipoallergenici per pelli sensibili_

—
—
—

Assortimento di 5 forme diverse:_

—
—

• 19 x 72 mm_

—
—

• 25 x 72 mm_

—
—

• 9 x 38 mm_

—
—

• 38 x 38 mm_

—
—
—

• Ø 12 mm_

—

Unità di vendita: 3200 (80 confezioni da 40 pezzi)