


**CATETERE FOLEY LATTICE 2 VIE ch/fr 16 - 30 ml -
sterile conf. 50 pz.**



CATETERI FOLEY STERILI IN LATTICE 100% RIVESTITI IN SILICONE - MONOUSO

- ch/fr 16

Cateteri a 2 vie - palloncino 30 ml - lunghezza 40 cm - valvola rigida

Scatola e manuale multilingue: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE,
FI, RO, BG, GR, HR, SI