


**CATETERE FOLEY LATTICE 3 VIE ch/fr 22 - 30 ml -
sterile conf. 50 pz.**



CATETERI FOLEY STERILI IN LATTICE 100% RIVESTITI IN SILICONE - MONOUSO

- ch/fr 22

Cateteri a 3 vie - palloncino 30 ml - lunghezza 40 cm - valvola rigida

Scatola e manuale multilingue: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE,
FI, RO, BG, GR, HR, SI