


**CATETERE FOLEY SILICONE 2 VIE ch/fr 18 - 30 ml -  
sterile conf. 10 pz.**



CATETERI FOLEY STERILI IN SILICONE 100% - MONOUSO - SENZA DEHP

Cateteri a 2 vie - palloncino 30 ml - lunghezza 40 cm - valvola rigida

- ch/fr 18

Scatola e manuale multilingue: GB, FR, IT, ES, PT, DE, NL, PL, CZ, SK, HU, EE, LT, LV, DK, SE,  
FI, RO, BG, GR, HR, SI