

STETOSCOPIO RIESTER DUPLEX DI ALLUMINIO - neonatale - blu



RIESTER DUPLEX® ALUMINIUM STETHOSCOPES

- Riester code: 4051
- DUPLEX® NEONATAL - blue
- Membrane: Ø 28 mm
- Length: 76 cm
- Weight: 106 g

Riester DUPLEX® stethoscopes combine first-class quality at an attractive price.

They offer outstanding acoustic performance.

Special membrane for perfect auscultation with particularly flat and soft non-chill rims for improved adaptation on the skin.

Bell Ø 36 mm.

Supplied with a pair of replacement eartips and a replacement membrane.

Made in Germany.