


CAVO ECG UNIVERSALE



CAVO ECG - 10 derivazioni Lunga durata, cavo testato singolarmente, compatibile con i seguenti ECG: BIONET: Cardiocare, Cardiotouch, Cardio 7 CARDIOLINE: E1, E1R, E3, Delta 300, Delta 360, C3. CARDIETTE (Elettronica Trentina): Autoruler 12/1, Autoruler 12/3, Excel 103, Excel 106, Astron Comp 1, Astron Comp 3, Astron Comp 6, Start 100 P, Start 100 H, Start 200, Start 200HV, AR600, AR 600 ADV, AR1200, AR1200ADV. CARDIORAPID CAMINA: K111, K112, K33, K36, K300 M, K300 P, K600 P, K600 D, K130 B, K 130 D. CONTEC: 100G, 300G, 600G, 1,200G EDAN: Smart SE-1, SE-3 ESAOTE: P80, P8000 GIMA: CARDIOGIMA 1M, 3M, 6M, 12, 12M HP: Pagewriter 1513, 154B/C, 1514D, 1515A/B, 1517A. KONTRON EKG P80 MACQUAIRE CT100, CT200. M4 MEDICAL: M traces NIHON - KOHDEN Cardiofax 651, 6581, 8110. PROGETTI Tutti i modelli: EPG junior, EPG view, EPG6 view, EPG6 view plus SCHILLER AT/1, AT/2, AT/3, AT/4, AT/5, AT/6, AT/10, CS6, CS12, CS100, CV6, P80.